

Light-Scape GT-1



Temperature Gradient Table
Very Bright RGB Light

Light-Scape GT-1


Introducing our newest instrument for scientific research

A temperature gradient table with adjustable high-power RGB lighting

This innovative table allows to maintain a temperature gradient, while easily changing light conditions to study microorganisms like algae, bacteria, seeds.

A simple yet powerful way to see how these organisms react to different environments, helping us learn more about the natural world.

- User defined light zones
- Adjustable height of the light fixture by stepper motors
- Solid design, easy to handle
- Dedicated PC software
- High brightness RGB light
- Adaptive active cooling of LEDs



Light-Scape GT-1 Technical specifications

Dimensios L x l x h (mm)	1300 x 1200 x 1200
Working area L x l (mm)	1200 x 600
Temperature gradient direction (cold - warm)	Left - Right along the 1200 dimension
Weight (kg)	200
Cold side temperature (°C)	-3
Warm side temperature (°C)	50
Power: 230 VAC, 50 Hz	Nominal 1500 W, Maximum 3300 W
Maximum white brightness (LUX)	50 000 under RGB strips, 100 000 under RGB power LEDs
Light zones - software defined	25 RGB addressable strips, 19 high power LEDs at the rear
LEDs height (mm)	25 - 200, adjustable from PC software

Light-Scape GT-1

Applications

The temperature and light gradient table offers a controlled environment for conducting experiments that require precise settings, making it an invaluable resource across a broad spectrum of disciplines.

1. Agricultural Research


- Seed Germination Studies: To understand the optimal conditions for seed germination and early growth stages, enabling the selection of the best planting times and conditions for various crops.
- Plant Growth and Development: Investigating how different light wavelengths and temperatures affect plant morphology, flowering time, and yield.

2. Microbiology

- Microorganism Growth Studies: Examining the growth rates of bacteria, fungi, and algae under varying temperatures and light conditions to understand their ecological roles and optimize conditions for beneficial strains.
- Antimicrobial Resistance Research: Testing how environmental stressors influence the development of resistance in microbes.

3. Environmental Science

- Climate Change Research: Simulating different environmental conditions to study plant and microorganism responses to global warming and altered light conditions due to atmospheric changes.



4. Pharmacology and Biochemistry

- Compound Stability Testing: Assessing how temperature and light affect the stability of pharmaceuticals and biochemical compounds.
- Enzyme Activity Studies: Investigating how varying environmental conditions influence the activity rates of enzymes.

5. Aquaculture

- Algae Cultivation: Optimizing conditions for growing algae used in biofuels, dietary supplements, and as a food source in aquaculture.

6. Educational Purposes

- Teaching Tool: Serving as a practical educational tool in classrooms and laboratories to teach students about the effects of environmental conditions on living organisms and chemical compounds.



CAMBUS LABS CUSTOM INSTRUMENTS

Manufacturing for industry since 1991

Contact us



+40736 34 34 54



www.cambus-labs.com



contact@cambus-labs.com

